


Building 103
Campbell Road
Stoke on Trent
ST4 4DE



Intelligent Orthopaedics

Decontamination

Product: Staffordshire Orthopaedic Reduction Machine (STORM)
Part Number: 1-0-0 (rev 1.0 onwards)

Details of product:

Re-useable device, used in the operating theatre to aid the reduction of lower limb fractures.

Sterility confirmation

Tested to HTM2010
Sterility confirmed using Porous-Load steriliser 134-137°C/3 minute cycle

Cleaning confirmation

Tested to HTM2030
Cleanliness confirmed using hand pre-wash and 6-minute automated cycle [HAMO type 1 (10 DIN) / type 2 (15 DIN)]

Approved by:

Date:

Dr C Ian Moorcroft
Technical Director
Intelligent Orthopaedics Ltd.

Intelligent Orthopaedics Ltd
Company Registration No: 4504094 Registered Office: Globe Court, Stoke on Trent, ST4 1ET