


Building 103
Campbell Road
Stoke on Trent
ST4 4DE



Intelligent Orthopaedics

Decontamination

Product: Staffordshire Orthopaedic Reduction Machine (STORM)
operating kit
Part Number: 1-10-0 (rev 1.0 onwards)

Details of product:

1 x 3.2 mm drill
2 x 4.5 mm cortical bone screw
2 x 2 mm tension wires

Sterility confirmation

Tested to HTM2010
Sterility confirmed using Porous-Load steriliser 134-137°C/3 minute cycle

Approved by:

Date:

Dr C Ian Moorcroft
Technical Director
Intelligent Orthopaedics Ltd.

Intelligent Orthopaedics Ltd
Company Registration No: 4504094 Registered Office: Glebe Court, Stoke on Trent, ST4 1ET